



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

	NHS-Biotin
Product Code	BB001
CAS No.	35013-72-0
Molecular Formula	$C_{14}H_{19}N_3O_5S$
Molecular Weight	341.38
Ex (nm)	----
Em (nm)	-----
Storage	4°C, protect from light

TEST	SPECIFICATION
Solubility:	DMSO
Quality:	HPLC(98+%)