



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

5-FITC

Product Code	DB001
CAS No.	3326-32-7
Molecular Formula	C ₂₁ H ₁₁ NO ₅ S
Molecular Weight	389.38
Ex (nm)	494
Em (nm)	520
Storage	Store in freezer (-20°C) and protect from light.

TEST

SPECIFICATION

Solubility:	DMSO
Quality:	HPLC(95+%)