



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

5-TRITC

Product Code	DB003
CAS No.	80724-19-2
Molecular Formula	C ₂₅ H ₂₁ N ₃ O ₃ S
Molecular Weight	443.52
Ex (nm)	543
Em (nm)	571
Storage	Store in freezer (-20°C) and protect from light.

TEST	SPECIFICATION
Solubility:	DMSO
Quality:	HPLC(95+%)