



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

6-TRITC

Product Code	DB004
CAS No.	80724-20-5
Molecular Formula	$C_{25}H_{21}N_3O_3S$
Molecular Weight	443.52
Ex (nm)	544
Em (nm)	572
Storage	Store in freezer (-20°C) and protect from light.

TEST

SPECIFICATION

Solubility:	DMSO
Quality:	HPLC(95+%)