



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

	DABSYL chloride
Product Code	DB005
CAS No.	56512-49-3
Molecular Formula	C ₁₄ H ₁₄ ClN ₃ O ₂ S
Molecular Weight	323.80
Ex (nm)	430
Em (nm)	----
Storage	Store in freezer (-20°C) and protect from light.

TEST	SPECIFICATION
------	---------------

Solubility:	Water
Quality:	HPLC(99+%)