



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

5-FAM

Product Code	DB007
CAS No.	76823-03-5
Molecular Formula	C ₂₁ H ₁₂ O ₇
Molecular Weight	376.32
Ex (nm)	492
Em (nm)	518
Storage	Store in freezer (-20°C) and protect from light.

TEST	SPECIFICATION
------	---------------

Solubility:	DMSO
Quality:	HPLC(97+%)