



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

---

## SPECIFICATION

6-FAM, SE

<b>Product Code</b>	DB011
<b>CAS No.</b>	92557-81-8
<b>Molecular Formula</b>	$C_{25}H_{15}NO_9$
<b>Molecular Weight</b>	473.39
<b>Ex (nm)</b>	495
<b>Em (nm)</b>	517
<b>Storage</b>	Store in freezer (-20°C) and protect from light.

### TEST

### SPECIFICATION

---

<b>Solubility:</b>	DMSO: $\geq 45$ mg/mL
<b>Quality:</b>	HPLC(95+%)