



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

---

## SPECIFICATION

5(6)-CR110,SE

<b>Product Code</b>	DB015
<b>CAS No.</b>	<b>254732-34-8</b>
<b>Molecular Formula</b>	$C_{25}H_{18}ClN_3O_7$
<b>Molecular Weight</b>	507.89
<b>Ex (nm)</b>	505
<b>Em (nm)</b>	525
<b>Storage</b>	Store in freezer (-20°C) and protect from light.

TEST	SPECIFICATION
------	---------------

---

<b>Solubility:</b>	DMSO
<b>Quality:</b>	HPLC(95+%)