



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

5(6)-TAMRA

Product Code	DB018
CAS No.	150347-56-1
Molecular Formula	$C_{25}H_{22}N_2O_5$
Molecular Weight	430.45
Ex (nm)	541
Em (nm)	565
Storage	Store in freezer (-20°C) and protect from light.

TEST	SPECIFICATION
------	---------------

Solubility:	DMSO
Quality:	HPLC(97+%)