



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

5(6)-TAMRA, SE

Product Code	DB021
CAS No.	246256-50-8
Molecular Formula	$C_{29}H_{25}N_3O_7$
Molecular Weight	527.53
Ex (nm)	546
Em (nm)	575
Storage	Store in freezer (-20°C) and protect from light.

TEST	SPECIFICATION
------	---------------

Solubility:	DMSO
Quality:	HPLC(95+%)