



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

5-TAMRA, SE

Product Code	DB022
CAS No.	150810-68-7
Molecular Formula	C ₂₉ H ₂₅ N ₃ O ₇
Molecular Weight	527.53
Ex (nm)	547
Em (nm)	575
Storage	Store in freezer (-20°C) and protect from light.

TEST	SPECIFICATION
Solubility:	DMSO
Quality:	HPLC(95+%)