



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

5(6)-ROX

Product Code	DB024
CAS No.	198978-94-8
Molecular Formula	$C_{33}H_{30}N_2O_5$
Molecular Weight	534.60
Ex (nm)	568
Em (nm)	591
Storage	Store in freezer (-20°C) and protect from light.

TEST	SPECIFICATION
------	---------------

Solubility:	DMSO
Quality:	HPLC(95+%)