



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

5(6)-ROX,SE

Product Code	DB027
CAS No.	114616-32-9
Molecular Formula	$C_{37}H_{33}N_3O_7$
Molecular Weight	631.70
Ex (nm)	576
Em (nm)	601
Storage	4°C, protect from light

TEST	SPECIFICATION
Solubility:	DMSO
Quality:	HPLC(95+%)