



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

---

## SPECIFICATION

5-ROX,SE

<b>Product Code</b>	DB028
CAS No.	209734-74-7
<b>Molecular Formula</b>	$C_{37}H_{33}N_3O_7$
<b>Molecular Weight</b>	631.70
<b>Ex (nm)</b>	573
<b>Em (nm)</b>	602
<b>Storage</b>	4°C, protect from light

TEST	SPECIFICATION
<b>Solubility:</b>	DMSO
<b>Quality:</b>	HPLC(95+%)