



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

AMCA,SE

Product Code	DB031
CAS No.	113721-87-2
Molecular Formula	$C_{16}H_{14}N_2O_6$
Molecular Weight	330.30
Ex (nm)	350
Em (nm)	448
Storage	4°C, protect from light

TEST

SPECIFICATION

Solubility:

DMSO

Quality:

HPLC(95+%)