



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

DABCYL acid

Product Code	DB038
CAS No.	6268-49-1
Molecular Formula	$C_{15}H_{15}N_3O_2$
Molecular Weight	269.30
Ex (nm)	-----
Em (nm)	-----
Storage	Store at room temperature

TEST	SPECIFICATION
------	---------------

Solubility:	DMF
Quality:	HPLC(97+%)