



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

DABCYL Acid, SE

Product Code	DB039
CAS No.	146998-31-4
Molecular Formula	$C_{19}H_{18}N_4O_4$
Molecular Weight	366.37
Ex (nm)	453
Em (nm)	-----
Storage	Store at room temperature

TEST

SPECIFICATION

Solubility:

DMF

Quality:

HPLC(98+%)