



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

APTS

Product Code	DB051
CAS No.	196504-57-1
Molecular Formula	$C_{16}H_8NNa_3O_9S_3$
Molecular Weight	523.40
Ex (nm)	424
Em (nm)	505
Storage	Store in freezer (-20°C) and protect from light.

TEST

SPECIFICATION

Solubility:	Water
Quality:	HPLC(98+%)