



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

Glycocholic acid sodium salt

Product Code	DE005
CAS No.	863-57-0
Molecular Formula	$C_{26}H_{42}NO_6Na$
Molecular Weight	487.61
CMC	13 mM(20-25°C)
Aggregation number	1.9
mol wt	micellar avg mol wt 900(lit.)
Storage	4°C, protect from light

TEST	SPECIFICATION
------	---------------

Quality:	HPLC(99+%)
-----------------	------------

Appearance:	White to off-white Solid
--------------------	--------------------------