



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

Sodium glycochenodeoxycholate

Product Code	DE008
CAS No.	16564-43-5
Molecular Formula	$C_{26}H_{42}NNaO_5$
Molecular Weight	471.61
CMC	----
Aggregation number	---
mol wt	-----
Storage	4°C, protect from light and moisture

TEST	SPECIFICATION
------	---------------

Quality:	HPLC(97+%)
-----------------	------------

Appearance:	White to off-white Solid
--------------------	--------------------------