



Hangzhou Handy Bio-technology Co., Ltd.

Tel: 0086-571-86691505

Fax: 0086-571-86691509

E-mail:sale@hzhandy.com

Website:www.hzhandy.com

SPECIFICATION

Isoluminol

Product Code	HG008
CAS No.	3682-14-2
Molecular Formula	$C_8H_7N_3O_2$
Molecular Weight	177.16
Ex (nm)	----
Em (nm)	-----
Storage	Store at 4°C. Moisture Sensitive.

TEST

SPECIFICATION

Solubility:	DMSO
Quality:	HPLC(98+%)